

ISO9001 & ISO13485 & ISO11137

Product Name:	Aseptic Connector				
Product Application:	NEST Aseptic Connector Collection provide a highly advantageous connection solution				
	for disposable bioprocess systems. Its unique design allows for aseptic connections, even				
	in non-sterile environments, ensuring both convenience and safety during fluid transfer				
	processes.				
	It is commonly used to connect the medium container to a bioreactor or storage bag				
	during cell culture, as well as for aseptic fluid transfer during feeding and sampling in				
	downstream processes of biopharmaceutical production.				
Product Materials:	PC (body) & medical-grade silicone (gasket)				
	Meets USP, Class 6 standards				
product Specifications:	Temperature range: stored at room temperature				
	Shelf life: 3 years after date of production (ensure package is in good)				
Sterilization:	No				
Pyrogens:	Non-Pyrogenic				
RNase/DNase:	DNase/RNase free				
Insoluble particles:	Reference "Chinese Pharmacopoeia", insoluble particle test method, ≤ 25 particles/mL for				
	particles ≥10µm, ≤3 particles/mL for particles ≥ 25µm.				
BSE/TSE/MGO:	No animal-derived or genetically modified ingredients or tissues, no TSE/BSE/GMO risks.				
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating				
0	procedures.				
	Melt film performance test: qualified Pressure resistance test: qualified				
	Radiation resistance test: qualified Freeze resistance test: qualified				
	Microbial challenge experiment: qualified Endotoxin test: qualified				
	Nuclease test: qualified Insoluble particulate test: qualified				
	High temperature and high pressure resistance test: qualified				
Features:	1. The gender-neutral design improves compatibility, making it easy to integrate into				
	disposable systems.				
	2. The intuitive three-step connection process (open, dock, pull) is simple and convenient,				
	eliminating the need for clamps, fixtures, or welding machines.				
	3. Supports E-beam sterilization (Max. 50 kGy) or autoclaving (121°C, 20 min).				
	4. The built-in silicone gasket withstands pressure up to 410 kPa at room temperature.				
	5. The main body is made of PC material, and the gasket is made of medical-grade silicone.				
	Both materials comply with USP Class VI standards.				
	6. Meets biocompatibility requirements of USP 87, USP 88, and ISO 10993 standards.				
	7. Undergoes physical tests and is suitable for demanding application environments.				
Precautions:	1. When connecting the two connectors, be careful not to flip the film. Keep the film in the				
	correct position as shown in the diagram, and then connect the two aseptic connectors.				
	2. When pulling and detaching the film, make sure to pull both pieces together. Do not pull				
	only one of them.				

DATA SHEET





ISO9001 & ISO13485 & ISO11137

Standard	Mini-type


Product Range

Aseptic Connector

Cat. No.	Product	Dia.	Packaging (Double Bagged)	
			/Pack	/Case
747521	Standard	1/4"(6.4 mm)	1	10
747531		3/8"(9.5 mm)	1	10
747541		1/2"(12.7 mm)	1	10
747551		3/4"(19 mm)	1	10
747611	Min-type	1/8"(3.2 mm)	1	10
747621		1/4"(6.4 mm)	1	10
747631		3/8"(9.5 mm)	1	10

Technical Drawing

Aseptic Connector:



Compatibility	Dimension				
Compatibility	A(mm)	B(mm)	C(mm)	D(mm)	
1/4"(standard)	45.37	36.2	16.3	8.1	
3/8" (standard)	45.37	37.8	17.75	10.94	
1/2" (standard)	45.37	41.45	22.25	14.2	
3/4" (standard)	45.37	51.8	33.1	20.5	
1/8" (min-type)	40.73	29.25	11.04	4.55	
1/4" (min-type)	40.73	33.44	15.23	8.1	
3/8" (min-type)	40.73	33.8	15.6	10.9	

 Head Office
 No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China

 Tel: +86+ 510-6800 6788
 Email: info@nest-wuxi.com

 Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix) NEST scientific 株式会社 (Yokohama,Japan) NEST Scientific Europe B.V (Netherlands) Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)