

ISO9001 & ISO13485 & ISO11137

Product Name:	PCR Scraper		
Product Application:	Used to compact the PCR sealing film, so that the PCR plate is sealed by the sealing		
	film to avoid liquid evaporation and contamination during the PCR amplification		
	process.		
Product Materials:	TPE		
product Specifications:	Temperature range: stored at room temperature		
	Shelf life: 5 Years after date of production (ensure package is in good)		
Sterilization:	No		
BSE/TSE:	These products are deemed animal free.		
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures.		
	Appearance inspection: pass Packaging inspection: pass		
Features	1. Made of TPE		
	2. Using with sealing film to make the sealing film fit better with the PCR plate.		

## **Product Range**

Cat.No.	Color	Spec (mm)	/Pack
411001	Blue	70*78	1

## **Product Picture**



Oversea

NEST USA (New Jersey/ Phoenix) NEST scientific 株式会社 (Yokohama,Japan) NEST Scientific Europe B.V (Netherlands) Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)