





ISO9001 & ISO13485 & ISO11137

Product Name:	PCR Sealing Film		
Product Application:	Nest PCR sealing film can be used with both Elisa plates and PCR plates. The film		
	can minimize the solution evaporation even under high temperature. Both common		
	and advanced films are available.		
Product Materials:	Polyolefin Film + pressure sensitive adhesive		
product Specifications:	Temperature range: stored at room temperature Shelf life: 3 Years after date of production (ensure package is in good)		
Sterilization:	No		
Pyrogens:	Non-Pyrogenic		
RNase/DNase:	DNase/RNase free		
BSE/TSE:	These products are deemed animal free.		
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating		
	procedures.		
	Appearance inspection: pass	Packaging inspection: pass	
Features	Common Sealing Film Product Features 1. Non-permeable film, medical grade glue		
	2. Compatible with most PCR plates		
	3. Good sealing performance		
	Advanced PCR Sealing Film Product Features		
	1. Non-tacky to skin and gloves		
	2. Minimal to no Auto-fluorescence		

Product Range

Cat. No.	Product Description	Spec (mm)	/Pack	/Case
410001	PCR Common Sealing Film.	146 × 81	100	500
410021	PCR Common Sealing Film, higher stickiness.	141 × 83	100	500
410011	Advanced sealing film for qPCR (high transparency).	141 × 78	100	500

DATA SHEET







ISO9001 & ISO13485 & ISO11137

Product Picture

