



Introduction


This product can be used in variety of lab oratory applications and work with a number of different sized tubes.

Feature




and easy operation.

Extremely compact



03





power supply.

Place the tube on

the top socket.

- Please connect the unit to the
- Press the tube slightly, and the



IP rating (DIN EN 60529)

Ambient Temp.

Power Supply

Power

Weight

Outside Dimensions

 Lift the tube, and the unit will stop shaking. Specifications	Mini Mixer NEST NEST
Cat. No.	105003
Oscillation Mode	Circumference Oscillation
Orbital Diameter	Ф4.5mm
Max. Capacity	80mL
Max. Tube Diameter	Ф30mm
Speed	3000rpm

IP40

5W

5~40 ℃

0.55kg

Ф96xH85mm

AC100~240V,50/60Hz