

ISO9001 & ISO13485 & ISO11137

Product Name:	Refill Pipette Tips					
Product Application:	1.Refill pipette tips, convenient for customers to use, reduce workload, shorten working					
	time, improve work efficiency.					
	2.Refill pipette tips are regular consumables for quantitative liquid transfer, used in the field					
	of biology, chemistry and other experimental direction with pipetting operation.					
Product Materials:	polypropylene (PP).					
	Meets USP, Class VI standards					
Product Specifications:	Temperature range: stored at room temperature					
	Shelf life: 5 Years after date of production (ensure package is in good)					
Sterilization:	Yes.					
	Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irrac					
	and dosimetrically released based on	ANSI/AAMI/ISO 11137				
Pyrogens:	Non-Pyrogenic					
RNase/DNase:	DNase/RNase free					
BSE/TSE:	These products are deemed animal free.					
Performance Test:	Each manufacturing lot is sampled and tested in accordance with standard operating					
	procedures.					
	Appearance inspection: qualified	Packaging inspection: qualified				
	Pipette fit test: qualified	Leak test: qualified				
Features	1.Both tips and boxes are autoclavable.					
	2.With precision mold and professional injection skill, Nest refill pipette tips has smooth					
	inner surface to ensure Low liquid attachment and accurate sample results during liquid					
	transfer.					
	3.Resistant to organic solvents.					
	4.ultra-clear polypropylene good for observation.					
	5.Nest refill pipette tips are compatible with most of the popular pipettes on the market.					
	6. Nest refill pipette tips is made of stacked brackets, different specifications of the suction					
	head inserted into the suction head box frame 5-10 stacked assembly, convenient for					
	customers to use, reduce costs, high utilization rate, only one box bottom, use up a box					
	frame and then take a box frame loading.					

Oversea

DATA SHEET





ISO9001 & ISO13485 & ISO11137


Product Range

Cat. No	Volume(µL)	Sterile	/Stack	/Pack	/Case
301081	10	No	96	960	9600
302081	200	No	96	960	9600
305081	300	No	96	480	4800
303081	1000	No	96	480	4800

Technical Drawing







 Head Office
 No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China

 Tel:
 +86+ 510-6800 6788
 Email: info@nest-wuxi.com


 Online:
 www.cell-nest.com
 Email: info@nest-wuxi.com


Oversea


NEST USA (New Jersey/ Phoenix) NEST scientific 株式会社 (Yokohama,Japan) NEST Scientific Europe B.V (Netherlands) Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)



ISO9001 & ISO13485 & ISO11137







 Head Office
 No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China

 Tel:
 +86+ 510-6800 6788
 Email: info@nest-wuxi.com

 Online:
 www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix) NEST scientific 株式会社 (Yokohama,Japan) NEST Scientific Europe B.V (Netherlands) Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)