






ISO9001 & ISO13485 & ISO11137

Product Name:	L-Spreader		
Product Application:	Ideal for spreading cells or bacteria evenly in the culture dish.		
Product Materials:	Polystyrene (PS)		
product Specifications:	: Temperature range: stored at room temperature  Shelf life: 5 Years after date of production (ensure package is in good		
Sterilization:	Yes.  Sterilized by E-beam, Sterility Assurance Level: SAL=10 <sup>-6</sup> . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137		
Pyrogens:	Non-Pyrogenic		
RNase/DNase:	DNase/RNase free		
BSE/TSE:	These products are deemed animal free.		
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating		
	procedures.		
	Appearance inspection: qualified Packaging inspection	n: qualified	
Features	Sterile L-Spreader, made of using medical PS.		
	<ol> <li>Smooth and flat surface that avoids problems such as contamination or uneven application of spread plate in bio-labs.</li> <li>Length of handle: 144mm, diagram of spreading: 33mm, suitable for 3.5cm, 6cm, 9d</li> <li>10cm petri dishes, etc. Ideal for spreading cells or bacteria evenly in the culture dish.</li> </ol>		

## **Product Range**

Cat. No.	Product Description	/Pack	/Case
711001	L-Spreader, Blue. Sterile	200	2000

## **Technical Drawing**



Head Office

No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China

Tel: +86+ 510-6800 6788 Email: <u>info@nest-wuxi.com</u>

Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix) NEST scientific 株式会社 (Yokohama, Japan) NEST Scientific Europe B.V (Netherlands) Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)