






ISO9001 & ISO13485 & ISO11137

Product Name:	Inoculating Loop & Needle				
Product Application:	Applied in microbiology experiments. For inoculating bacteria on agar medium.				
Product Materials:	Polypropylene (PP)				
product Specifications:	Temperature range: stored at room temperature				
•	Shelf life: 5 Years after date of production (ensure package is in good)				
Sterilization:	Yes.				
Pyrogens:	Non-Pyrogenic				
RNase/DNase	DNase/RNase free				
BSE/TSE:	These products are deemed animal free.				
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures.				
	Appearance inspection: qualified Packaging inspection: qualified				
Features	1. PP material, hydrophilic.				
	2. Different colors for identification.				
	3. Flexible handles for easier collection and inoculation.				
	5. Easy-tear sterile individual packaging.				

Product Range

Cat. No.	Volume(μL)	Name	Color	/Pack	/Case
716001	1	Inoculating Needle	White	400	4000
717101	1	Inoculating Loop	Blue	400	4000
718202	10	Inoculating Loop	Yellow	400	4000

Product Picture










ISO9001 & ISO13485 & ISO11137

Technical Drawing







Inoculating Needle

Inoculating Loop(1 µL)

Inoculating Loop(10 µL)