

<b>Product Name:</b>	<b>Cryo Boxe</b>
<b>Product Application:</b>	<b>For centralized placement of cryogenic vials.</b>
<b>Product Materials:</b>	<b>PC (cap) + ABS (bottom)</b>
<b>product Specifications:</b>	Temperature range: stored at room temperature  Shelf life: 5 Years after date of production ( ensure package is in good )
<b>Sterilization:</b>	No
<b>BSE/TSE:</b>	These products are deemed animal free.
<b>Performance Testing:</b>	Each manufacturing lot is sampled and tested in accordance with standard operating procedures.  Appearance inspection: pass      Packaging inspection: pass
<b>Features</b>	<ol style="list-style-type: none"> <li>1. The cap is made of polycarbonate plastic (PC), and the bottom is made of acrylonitrile butadiene styrene plastic (ABS).</li> <li>2. Clear polycarbonate cover for easy identification.</li> <li>3. Individually packed.</li> <li>4. Alphanumeric grid for easy sample identification.</li> <li>5. Each box includes drain holes and air vents to minimize condensation.</li> <li>6. Withstand temperature range: Approx.-196°C to 121°C.</li> </ol>


**Product Range**

Cat. No.	Product Description	/Pack	/Case
616224	For 0.5 ml vials	9*9	56
616325	For 1.2 ml vials	9*9	48
616021	For 2.0 ml vials	9*9	40
616421	For 3.0 ml vials	9*9	32
616121	For 4.0 ml vials	9*9	32
616521	For 5.0 ml vials	9*9	32
616204	For 0.5ml vials, Internal Thread Cryogenic Vial Only	10*10	56
616305	For 1.2ml vials, Internal Thread Cryogenic Vial Only	10*10	48
616001	For 2.0ml vials, Internal Thread Cryogenic Vial Only	10*10	40


**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
 Tel: +86+ 510-6800 6788      Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
 Online: [www.cell-nest.com](http://www.cell-nest.com)

**Oversea** NEST USA (New Jersey/ Phoenix)  
 NEST scientific 株式会社 (Yokohama, Japan)  
 NEST Scientific Europe B.V (Netherlands)  
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

Product Picture







Technical Drawing



**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
 Tel: +86+ 510-6800 6788 Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
 Online: [www.cell-nest.com](http://www.cell-nest.com)

**Oversea** NEST USA (New Jersey/ Phoenix)  
 NEST scientific 株式会社 (Yokohama, Japan)  
 NEST Scientific Europe B.V (Netherlands)  
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

ISO9001 & ISO13485 & ISO11137




**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
 Tel: +86+ 510-6800 6788 Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
 Online: [www.cell-nest.com](http://www.cell-nest.com)


**Oversea** NEST USA (New Jersey/ Phoenix)  
 NEST scientific 株式会社 (Yokohama, Japan)  
 NEST Scientific Europe B.V (Netherlands)  
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

ISO9001 & ISO13485 & ISO11137


**10\*10 Cryo Boxes for 0.5 ml Vials**  
Internal Thread Cryogenic Vial Only



**10\*10 Cryo Boxes for 1.2 ml Vials**  
Internal Thread Cryogenic Vial Only



**10\*10 Cryo Boxes for 2.0 ml Vials**  
Internal Thread Cryogenic Vial Only



**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
Tel: +86+ 510-6800 6788 Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
Online: [www.cell-nest.com](http://www.cell-nest.com)

**Oversea** NEST USA (New Jersey/ Phoenix)  
NEST scientific 株式会社 (Yokohama, Japan)  
NEST Scientific Europe B.V (Netherlands)  
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)