

Product Name:	Nitrile Powder Free Examination Gloves		
Product description:	Nitrile Gloves, Powder Free Non sterile (Non Medical/Industrial)		
Product Materials:	Synthetic Nitrile Butadiene Rubber (NBR)		
Feature:	Powder free, online chlorinated		
Design:	Ambidextrous, standard beaded cuff and fingertip		
Color	Black		
Product packaging mode:	100 gloves per box by weight and 10 boxes per carton <i>*Subject to customer requirement</i>		
Physical Dimension:	MSH	ASTM D6319	
Glove length (mm)	≥ 240	S -220, min, M, L, XL -230, min	
Glove Weight (g)	MSH	ASTM D6319	
Small	6.2 ± 0.2	NA	
Medium	6.5 ± 0.2		
Large	6.8 ± 0.2		
Extra large	7.1 ± 0.2		
Width (mm)	MSH	ASTM D6319	
Small	80 ±10	80 ±10	
Medium	95 ±10	95 ±10	
Large	110 ±10	110 ±10	
Extra large	120 ±10	120 ±10	
Glove's thickness (mm), (measured in single wall):			
	MSH	ASTM D6319	
Finger (13 ±3mm from fingertip)	0.15, min (5.9 mil)	0.05, min	
Palm (Approximate center of palm, with respect to glove size)	0.09, min (3.5 mil)	0.05, min	
Cuff (Approximate 25±5mm from the cuff edge)	0.08, min (3.1 mil)	NA, min	
Physical Properties:	<p>Tensile strength (MPa):</p> <ul style="list-style-type: none"> • Before Aging: Min 14 • After challenge testing (Aging 7 days at 70 ± 2°C): Min 14 <p>Elongation (%):</p> <ul style="list-style-type: none"> • Before Aging: Min 500 • After challenge testing (Aging 7 days at 70 ± 2°C): Min 400 <p>Force at Break (N):</p> <ul style="list-style-type: none"> • Before Aging: NA • After challenge testing (Aging 7 days at 70 ± 2°C): NA 		

Product Performance Requirements:	Characteristic	Inspection level	AQL
	1000mL Watertightness	GI	4.0
	Physical dimensions	S2	4.0
	Physical properties	S2	4.0
	Visual defects - Critical	GI	0.65
	Visual defects – Major & Minor	GI	6.5

Product Range

Cat.No.	Size	Product	/Pack	/Case
902801	Extra Small	Extra Thick Nitrile examination gloves, Black, 5 mil, powder free, Non-Sterile,	100	1,000
902811	Small	Extra Thick Nitrile examination gloves, Black, 5 mil, powder free, Non-Sterile,		
902821	Medium	Extra Thick Nitrile examination gloves, Black, 5 mil, powder free, Non-Sterile,		
902831	Large	Extra Thick Nitrile examination gloves, Black, 5 mil, powder free, Non-Sterile,		
902841	Extra Large	Extra Thick Nitrile examination gloves, Black, 5 mil, powder free, Non-Sterile,	90	900

Picture



Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
 Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
 Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
 NEST scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)