

Pipette Brand	Series	Channel	Volume (µL)	301006	305006/ 302106	304006/ 303206	301016/ 301081	305016/302116/ 302081/305081	304016/ 303216/30308	
Axygen	Axyper	Single-channel	0.1-2	•			•			
			0.5-10	•			•			
			2-20		•			•		
			5-50		•			•		
			10-100		•			•		
			20-200		•			•		
		100-1000			•			•		
		Multi-channel	1-10	•				•		
			5-50		•			•		
			20-200		•			•		
50-300			•			•				
Biohit	Proline	Single-channel	0.1-2.5	•			•			
			0.2-10	•			•			
			0.5-10	•			•			
			5-50		•			•		
			5-100		•			•		
			10-250		•			•		
			10-500			•			•	
			20-200		•			•		
			50-1000			•			•	
			100-1000			•			•	
	Multi-channel	0.2-10	•				•			
		5-100		•			•			
		25-250		•			•			
		50-250		•			•			
		5-100		•			•			
	e-Line	Single-channel	0.1-10	•			•			
			0.2-10	•			•			
			5-120		•			•		
			10-300		•			•		
		50-1000			•			•		
Multi-channel		0.2-10	•				•			
		5-120		•			•			
		10-300		•			•			
m-Line	Single-channel	0.1-3	•			•				
		0.5-10	•			•				
		10-100		•			•			
		20-100		•			•			
		100-1000			•			•		
	Multi-channel	0.5-10	•				•			
		5-100		•			•			
BrandTech	Transferpete	Single-channel	0.1-1	•			•			
			0.5-10	•			•			
			2-20		•			•		
			5-50		•			•		
			10-100		•			•		
			20-200		•			•		
			25-250		•			•		
			100-1000			•			•	
		Multi-channel	0.5-10	•				•		
			2-20		•			•		

		Multi-channel	2.5-25		•		•	
			5-50		•		•	
			10-100		•		•	
			20-200		•		•	
			30-300		•		•	
Capp Aero		Single/Multi-channel	0.2-2	•		•		
			0.5-10	•		•		
			5-50		•		•	
			10-100		•		•	
			25-200		•		•	
			50-300		•		•	
100-1000			•		•			
Labnet	Labpette Discovery or Biopette™ Plus	Single-channel	0.1-2	•		•		
			0.5-10	•		•		
			2-20		•		•	
			10-100		•		•	
			20-200		•		•	
			100-1000			•		•
	Multi-channel	1-10	•		•			
		5-50		•		•		
		20-200		•		•		
		50-300		•		•		
	BioPette E	Single-channel	0.5-10	•		•		
			2-20	•		•		
			10-200		•		•	
100-1000					•		•	
Multi-channel		0.5-10	•		•			
		2-20	•		•			
	10-200		•		•			
100-1200			•		•			
Matrix	Impact/Impact ²	Single/Multi-channel	2-125		•		•	
Nichryo/Oxford	Benchmate	Single-channel	0.5-10	•		•		
			2-20		•		•	
			10-50		•		•	
Rainin	EDP3	Single-channel	E-MIC-10,0.5-10	•		•		
			E2-25,2.5-25		•		•	
			E2-100,10-100		•		•	
			E2-250,25-250		•		•	
			E2-1000,100-1000			•		•
	Multi-channel	EP-M8-10,0.5-10	•		•			
		EP-M8-250,10-250		•		•		
	Latch Mode	Single-channel	R2,0.1-2	•		•		
			R10,0.2-10	•		•		
			R20,2-20		•		•	
R100,10-100				•		•		
R200,20-200				•		•		
R1000,100-1000					•		•	
Socorex	Callbra/Accura	Single-channel	0.1-2	•		•		
			0.2-2	•		•		
			0.5-10	•		•		
			1-10	•		•		
			2-20		•		•	
			10-100		•		•	
			20-200		•		•	
			100-1000			•		•

		Multi-channel	1-10	•			•			
			5-50						•	
			10-100						•	
			20-200						•	
			50-200						•	
VWR	VWR and VWR Signature	Single-channel	0.1-2	•			•			
			0.5-10	•			•			
			2-20						•	
			5-50						•	
			10-100						•	
			20-200						•	
			50-250						•	
		100-1000				•			•	
		Multi-channel	1-10	•				•		
			5-50							•
20-200								•		
50-300								•		
Eppendorf	Eletronic	Single/Multi-channel	0.5-10	•			•			
			5-100						•	
			2-200						•	
			50-1000					•		•
	Reference	Single/Multi-channel	0.1-2.5	•			•			
			0.1-10	•			•			
			0.5-10	•			•			
			2-20						•	
			10-100						•	
			20-200						•	
			50-200						•	
	100-1000					•		•		
	Research	Single-channel	0.1-2.5	•			•			
			0.5-10	•			•			
			2-20						•	
			10-100						•	
			20-200						•	
		100-1000					•		•	
		Multi-channel	0.5-10	•				•		
			5-100							•
			10-100							•
			50-100							•
	20-300								•	
	Research Pro	Single-channel	20-300						•	
			50-1000					•		•
		Multi-channel	0.5-10	•				•		
	20-300		•				•			
	Research Plus	Single/Multi-channel	100-1000				•			•
			Single-channel	0.2-2	•			•		
				0.3-3	•			•		
0.2-10				•			•			
0.5-5				•			•			
0.5-10				•			•			
1-10				•			•			
5-40									•	
5-50									•	

Thermo Scientific	Finnpipette		10-100		•			•	
			40-200		•			•	
			30-300		•			•	
			100-1000			•			•
			200-1000			•			•
	Multi-channel	1-10	•			•			
		5-50		•			•		
		50-300		•			•		
	Finnpipette Digital	Single-channel	40-200		•				
	Finnpipette electronic	Multi-channel	1-10	•			•		
30-300				•			•		
100-1000					•			•	
Gilson Pipetman	Pipetman Classic	Single-channel	P2,0,1-2	•			•		
			P2,0,2-2	•			•		
			P10,0,2-10	•			•		
			P20,2-20		•			•	
			P100,10-100		•			•	
			P200,20-200		•			•	
			P1000,100-1000			•			•
	Pipetman Neo	Single-channel	P2N,0,2-2	•			•		
			P10N,1-10	•			•		
			P20N,2-20		•			•	
			P100N,10-100		•			•	
			P100N,20-100		•			•	
		P200N,20-200		•			•		
	Multi-channel	P8×20N,2-20		•			•		
	Pipetman Ultra	Single-channel	U200,20-200		•			•	
		Multi-channel	P200-M8,20-200		•			•	



Pipette Tips Compatibility Reference Tool (with Filter)

Pipette Brand	Series	Channel	Volume (µL)	311001	312001/ 316001/310012	313001	311012	316012/ 312012	313012	314016	310012	313006/31 3016		
Axygen	AxyPET	Single-channel	0.1-2	•			•			•				
			0.5-10	•			•			•				
			2-20		•			•				•		
			5-50		•			•						
			10-100		•			•						
			20-200		•			•						
		100-1000				•			•			•		
		Multi-channel	1-10	•				•			•			
			5-50		•			•						
			20-200		•			•						
50-300			•			•								
Biohit	Proline	Single-channel	0.1-2.5	•			•			•				
			0.2-10	•			•			•				
			0.5-10	•			•			•				
			5-50		•			•						
			5-100		•			•						
			10-250		•			•						
			10-500				•			•				
			20-200		•			•						
			50-1000				•			•				
			100-1000				•			•			•	
	Multi-channel	0.2-10	•				•			•				
		5-100		•			•							
		25-250		•			•							
		50-250		•			•							
		5-100		•			•							
		5-100		•			•							
	e-Line	Single-channel	0.1-10	•			•				•			
			0.2-10	•			•				•			
			5-120		•			•						
			10-300		•			•						
50-1000						•			•			•		
Multi-channel		0.2-10	•				•			•				
		5-120		•			•							
		10-300		•			•							
		5-100		•			•							
		5-100		•			•							
m-Line	Single-channel	0.1-3	•			•				•				
		0.5-10	•			•				•				
		10-100		•			•							
		20-100		•			•							
	Multi-channel	100-1000				•			•			•		
		0.5-10	•				•			•				
		5-100		•			•							
		30-300		•			•							
BrandTech	Transferpet te	Single-channel	0.1-1	•			•			•				
			0.5-10	•			•			•				
			2-20		•			•			•			
			5-50		•			•						
			10-100		•			•						
			20-200		•			•						
			25-250		•			•						
			100-1000				•			•			•	
			Multi-channel	0.5-10	•				•			•		
				2-20		•			•				•	
	2.5-25			•			•							
	5-50			•			•							
	Single/Multi -channel	10-100		•			•							
		20-200		•			•							
30-300			•			•								
0.2-2		•				•			•					
0.5-10		•				•			•					
5-50			•			•								
Capp Aero	Single/Multi -channel	10-100		•			•							
		25-200		•			•							

			50-300			•			•						
			100-1000				•			•			•		
Labnet	Labpette Discovery or Biopette™ Plus	Single-channel	0.1-2	•				•			•				
			0.5-10	•				•			•				
			2-20			•			•				•		
			10-100			•			•						
			20-200			•			•						
			100-1000					•			•			•	
	Multi-channel	1-10	•					•			•				
		5-50			•				•						
		20-200			•				•						
		50-300			•				•						
		BioPette E	Single-channel	0.5-10	•				•			•			
		2-20		•					•			•	•		
	10-200				•				•						
	Multi-channel	100-1000					•				•			•	
0.5-10		•						•			•				
2-20		•						•			•	•			
10-200				•				•							
Matrix	Impact/Impact ²	Single/Multi-channel	2-125			•			•						
			Nichryo/Oxford	Benchmate	Single-channel	0.5-10	•				•			•	
			2-20					•			•			•	
			10-50					•			•				
Rainin	EDP3	Single-channel	E-MIC-10,0.5-10	•					•				•		
			E2-25,2.5-25			•				•					
			E2-100,10-100			•				•					
			E2-250,25-250			•				•					
			E2-1000,100-1000					•				•		•	
	Latch Mode	Single-channel	EP-M8-10,0.5-10	•					•				•		
			EP-M8-250,10-250			•				•					
			R2,0.1-2	•						•				•	
			R10,0.2-10	•						•				•	
			R20,2-20			•				•				•	
Socorex	Callbra/Accura	Single-channel	R100,10-100			•			•						
			R200,20-200			•			•						
			R1000,100-1000					•				•		•	
			0.1-2	•					•			•			
			0.2-2	•					•			•			
VWR		VWR and VWR Signature	Single-channel	0.5-10	•					•			•		
				2-20			•				•			•	
				5-50			•				•				
				10-100			•				•				
				20-200			•				•				
	Multi-channel		50-250			•				•					
			100-1000					•				•		•	
			1-10	•						•			•		
			5-50			•				•					
			20-200			•				•					
Eletronic	Single/Multi-channel	50-300			•				•						
		0.5-10	•						•			•			
		5-100			•				•						
			2-200			•			•						

Eppendorf	Reference	Single/Multi-channel	50-1000			•			•			•		
			0.1-2.5	•			•			•				
			0.1-10	•			•			•				
			0.5-10	•			•			•				
			2-20			•			•			•		
			10-100			•			•					
			20-200			•			•					
			50-200			•			•					
	100-1000					•			•			•		
	Research	Single-channel	0.1-2.5	•			•			•				
			0.5-10	•			•			•				
			2-20			•			•			•		
			10-100			•			•					
			20-200			•			•					
		100-1000					•			•			•	
		Multi-channel	0.5-10	•				•			•			
			5-100			•				•				
			10-100			•				•				
			50-100			•				•				
			20-300			•				•				
30-300					•				•					
50-1200						•			•			•		
Research Pro	Single-channel	20-300			•			•						
		50-1000					•		•					
	Multi-channel	0.5-10	•			•			•					
		20-300	•			•			•					
Research Plus	Single/Multi-channel	100-1000				•			•			•		
Thermo Scientific	Finnpipette	Single-channel	0.2-2	•			•			•				
			0.3-3	•			•			•				
			0.2-10	•			•			•				
			0.5-5	•			•			•				
			0.5-10	•			•			•				
			1-10	•			•			•				
			5-40			•			•					
			5-50			•			•					
			10-100			•			•					
			40-200			•			•					
			30-300			•			•					
			100-1000					•			•			•
	200-1000					•			•			•		
	Finnpipette Digital	Single-channel	1-10	•			•			•				
			5-50			•			•					
50-300					•			•						
Finnpipette electronic	Multi-channel	1-10	•			•			•					
		30-300			•			•						
		100-1000					•			•			•	
Gilson Pipetman	Pipetman Classic	Single-channel	P2,0.1-2	•			•			•				
			P2,0.2-2	•			•			•				
			P10,0.2-10	•			•			•				
			P20,2-20			•			•			•		
			P100,10-100			•			•					
			P200,20-200			•			•					
	Pipetman Neo	Single-channel	P1000,100-1000				•			•			•	
			P2N,0.2-2	•			•			•				
			P10N,1-10	•			•			•				
			P20N,2-20			•			•			•		
			P100N,10-100			•			•					
			P100N,20-100			•			•					
	Pipetman Ultra	Multi-channel	P200N,20-200			•			•					
			P8x20N,2-20			•			•					
			U200,20-200			•			•					
		P200-M8,20-200			•			•						