

<b>Product Name:</b>	<b>Cell Scrapers</b>
<b>Product Application:</b>	<b>Specially designed Nest® cell scrape can be used to scrap and collect cells easily and effectively. Push the handle with a slight pressure to change the blade angle. Then push down the handle towards the bottom of the container. Rotate the handle slightly to make the blade twist to right direction.</b>
<b>Product Materials:</b>	<b>HIPS (Blade) + PP (handle)</b>
<b>product Specifications:</b>	Temperature range: stored at room temperature Shelf life: 5 Years after date of production (ensure package is in good)
<b>Sterilization:</b>	<b>Yes.</b> Sterilized by E-beam, Sterility Assurance Level: SAL=10 <sup>-6</sup> . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
<b>Pyrogens:</b>	Non-Pyrogenic
<b>RNase/DNase:</b>	DNase/RNase free
<b>BSE/TSE:</b>	These products are deemed animal free.
<b>Performance Testing:</b>	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified      Packaging inspection: qualified
<b>Features</b>	<ol style="list-style-type: none"> <li>1. Blade is made of HIPS and handle is made of PP</li> <li>2. Free rotating blade, make the blade twist to right direction</li> <li>3. Easy-tear sterile packing</li> <li>4. Sterilized by E-beam, SAL=10<sup>-6</sup></li> <li>5. Non-pyrogenic and DNase/RNase free</li> </ol>



**Product Range**

Cat. No.	Spec(mm)		Application	/Case
	Total Length	Blade Width		
710001	220	13	blue handle & transparent blade	100
710011	280	20		

**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
Tel: +86+ 510-6800 6788      Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
Online: [www.cell-nest.com](http://www.cell-nest.com)

**Oversea** NEST USA (New Jersey/ Phoenix)  
NEST scientific 株式会社 (Yokohama, Japan)  
NEST Scientific Europe B.V (Netherlands)  
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

Technical Drawing



**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
Tel: +86+ 510-6800 6788 Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
Online: [www.cell-nest.com](http://www.cell-nest.com)

**Oversea** NEST USA (New Jersey/ Phoenix)  
NEST scientific 株式会社 (Yokohama, Japan)  
NEST Scientific Europe B.V (Netherlands)  
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)