





ISO9001 & ISO13485 & ISO11137

Product Name:	PCR Sealing Film				
Product Application:	Nest PCR sealing film can be used with both Elisa plates and PCR plates. The film				
	can minimize the solution evaporation even under high temperature. Both common				
	and advanced films are available.				
Product Materials:	Polyolefin Film + pressure sensitive adhesive				
product Specifications:	Temperature range: stored at room temperature				
	Shelf life: 5 Years after date of production (ensure package is in good)				
Sterilization:	No				
Pyrogens:	Non-Pyrogenic				
RNase/DNase:	DNase/RNase free				
BSE/TSE:	These products are deemed animal free.				
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating				
	procedures.				
	Appearance inspection: pass	Packaging inspection: pass			
Features	Common Sealing Film Product Features				
	1. Non-permeable film, medical grade	glue			
	2. Compatible with most PCR plates				
	3. Good sealing performance				
	Advanced PCR Sealing Film Product Features				
	1. Non-tacky to skin and gloves				
	2. Minimal to no Auto-fluorescence				

Product Range

Cat. No.	Product Description	Spec (mm)	Usable Film Dimension	/Pack	/Case
410001	PCR Common Sealing Film.	146 × 81	124×81	100	500
410021	PCR Common Sealing Film, higher stickiness.	141 × 83	118×80	100	500
410011	Advanced sealing film for qPCR (high transparency).	141 × 78	116×78	100	500

DATA SHEET







ISO9001 & ISO13485 & ISO11137

Product Picture

