

ISO9001 & ISO13485 & ISO11137

Product Name:	Cover Glasses			
Product Application:	Suitable for a variety of cell staining, including immunohistochemistry and immunofluorescence staining followed by laser confocal observation.			
Product Materials:	High quality borosilicate glass			
product Specifications:	Temperature range: stored at room temperature			
	Shelf life: 3 Years after date of production ( ensure package is in good )			
Sterilization:	Yes.			
Pyrogens:	Non-pyrogenic			
RNase/DNase:	DNase/RNase free			
BSE/TSE:	These products are deemed animal free.			
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating			
	procedures.			
	Appearance inspection: qualified Cell culture test: qualified			
	Packaging inspection: qualified			
Features	1. High quality borosilicate glass, resistant to organic solvents.			
	2. No pyrogen, no cytotoxicity, DNase/RNase free			
	3. Autoclavable			
	4. Not TC treated.			
	5. Thickness: 0.16-0.19mm			
	6. Can be used in confocal experiment.			
	7. Each package has an independent batch number for easy quality tracking and			
	traceability.			
Packages	Each package has a spare quantity to balance the loss during transportation and use.			

## **Product Range**

Cat. No.	Glass Diameter(mm)	Thickness (mm)	Application	/Pack
801010	14	0.16±0.02	24 Well Cell Culture Plate	100
801007	15	0.16±0.02	24 Well Cell Culture Plate	100
801011	18	0.16±0.02	12 Well Cell Culture Plate	100
801008	20	0.16±0.02	12 Well Cell Culture Plate	100
801009	25	0.16±0.02	6 Well Cell Culture Plate	100

 Head Office
 No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China

 Tel: +86+ 510-6800 6788
 Email: info@nest-wuxi.com

 Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix) NEST scientific 株式会社 (Yokohama,Japan) NEST Scientific Europe B.V (Netherlands) Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

## **DATA SHEET**



## **Product Picture**





No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Head Office Jiangsu, China Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com Online: www.cell-nest.com

NEST USA (New Jersey/ Phoenix) Oversea

NEST scientific 株式会社 (Yokohama,Japan) NEST Scientific Europe B.V (Netherlands) Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)