

Product Name:	NEST Automation Tips		
Product Application:	It can be used for research and development and other commonly used		
	applications with needs of high-throughput pipetting in genomics, proteomics,		
	cytomics, immunoassays, metabolomics, biopharmaceuticals.		
	1. NEST Automatic Tips have two types: clear and black.		
	2. Tecan automation tips may be used with the TECAN Freedom EVOlyzers and		
	fully automatic sampling systems.		
	3. Hamilton automation conductive tips are used in conjunction with Hamilton		
	series workstations.		
	4. Beckman robotic tips are used in Beckman FX/NX,Biomek i5/ i7 and Biomek		
	3000.		
Product Materials:	Medical-grade polypropylene (PP)		
product Specifications:	Temperature range: stored at room temperature		
	Shelf life: 5 Years after date of production (ensure package is in good)		
Sterilization:	Yes.		
	Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated		
	and dosimetrically released based on ANSI/AAMI/ISO 11137.		
Pyrogens:	Non-Pyrogenic		
RNase/DNase:	DNase/RNase free		
BSE/TSE:	These products are deemed animal free.		
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating		
	procedures.		
	Appearance inspection: qualified Air tightness test: qualified		
	Tests for resistance and CV values: qualified Packaging inspection: qualified		
Features	1. Medical-grade polypropylene (PP) materials is Imported and selected strictly		
	accordance with the industry requirements to ensure the stability of materials between		
	batches.		
	2. The production is produced automatically in one hundred thousand level purification		
	workshops, to ensure no DNA, RNA, proteases and pyrogens. Tips are packed in separate		
	cavities with uniform conductivity: tips of each box share the same cavity number to ensure		
	the product traceability and the conductivity uniformity and to improve the accuracy of		
	experimental results, greatly reducing the difference between individual products.		
	3. The process technology is unique, to ensure that the Tips has a smooth inner surface,		
	greatly reducing the amount of liquid residue.		
	4. Compatible with a variety of automated sample processing platforms. The structure		
	drawing test is according to the original adaptor, with a mature injection molding process, to		
	ensure that the product has good air tightness and adaptability, improve the product in the		
	work of mechanical accuracy.		
	 The control is strictly, to according to customer's needs for effective quality inspection, 		
	with no DNA, RNA, protease and heat source testing, sterility testing, etc.		
	אונד הי בוית, דוות, דוסנכמש מוע ווכמו שטווכב ובשנווע, שלווונע ובשנווע, בנט.		

 Head Office
 No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China

 Tel:
 +86+ 510-6800 6788
 Email:
 info@nest-wuxi.com

 Online:
 www.cell-nest.com
 info@nest-wuxi.com

DATA SHEET



	ISO9001 & ISO13485 & ISO11137
	6. Air tightness detection: the precision detector is used to check the air tightness
	between the Tips and the adapter, to ensure that each batch of products have good air
	tightness.
	7. The test equipment is used online to control the size and concentricity of products in
	production line, saving time and improving product quality.
	8. There are two styles with and without filter. The Automatic Filter Tips are filtered
	through porous tissue to ensure optimal performance and is extremely hydrophobic,
	forming a robust barrier against aerosols and eliminating the risk of cross-contamination of
	the sample.
	9. Strong package: the high-strength blister card outer package with thickened wall,
	which is impact-resistant and dropping-resistant, ensures the integrity and safety of the
	product under harsh transport conditions.
	NEST Automatic Tips have two types
	NEST Automation Conductive Tips:
	NEST Automation black Conductive Tips are made of Medical-grade black conductive
	polypropylene (PP) materials
	· Tests for resistance and CV values: a certain number of products are sampled for each
	batch for inspection, which is conducted using precise and unique measuring tools to
	ensure the accuracy and reliability of the test results to ensure the uniform conductivity
	of products.
	NEST Automatic Transparent Tips:
	· NEST Automatic transparent Tips are made of polymer material polypropylene (PP),
	and the surface is treated with a special process to have high hydrophobicity.
	· The automatic transparent tip is carefully made for special purposes and is designed
	to bring more efficient efficiency and repeatability to your work.
Types of Packaging	There are two packaging methods: Blister Package and Plastic Rack
	Blister Package:
	Sheet is a new kind of materials developed in recent years, with its good toughness, high
	strength, high transparency, and can be recycled, combustion does not produce harmful
	gases, excellent environmental performance and convenient user use, suction head new
	NEST conductive immovable blister, blister process will be transparent plastic plate made
	specific convex shape of transparent plastic, cover on the surface of the product, protect,
	and beautify the role of the product.
	1. It has good transparency and smoothness, and good display effect.
	2. With good mechanical strength, bubble shell material outer packaging, wall thickness
	size, anti-impact, anti-drop, to ensure the integrity and safety of the product under harsh
	transportation conditions.
	3. Good barrier performance for oxygen and water vapor
	4. Good chemical resistance can withstand a variety of chemical material erosion.
	5. Non-toxic, reliable health performance can be used for food, medicine, and medical

DATA SHEET



ISO9001 & ISO13485 & ISO11137

equipment packaging, and can tolerate gamma ray disinfection of its packaging items.

6. Good adaptability to environmental protection, can be economical and convenient recycling, do not produce harmful substances harmful to the environment.

7. The thermal forming performance is good, and the ordinary sheet material, in line with the environmental protection requirements of an ideal packaging material.

8. Convenient for users in the automation of large quantities of high efficiency.

Plastic Racked:

1. The plastic box is made of high-strength PP material, which has good strength and is not easy to deform.

2. The product size conforms to the industry SBS size standard and is suitable for a variety of automated pipetting workstation platforms.

3. Plastic boxed suction heads solve the unmanned contact needs of special application

enterprises, avoid hand contact, reduce pollution risks, and improve work efficiency.

4. No DNA\RNase, no endotoxin

Product Range

Automation Conductive Tips for Tecan, Black			
Images	Cat. No.	Product Description	Packing
		50 µl Automation Tips for Tecan,	96 Tips*2 Blisters/pk, 2304/cs
	332002	Conductive, Black, Sterile, Blister	
		Packed	
		50 µl Automation Filter Tips for Tecan,	96 Tips*2 Blisters/pk, 2304/cs
	332012	Conductive, Black, Sterile, Blister	
		Packed,	
		200 µl Automation Tips for Tecan,	96 Tips*2 Blisters/pk, 2304/cs
	332102	Conductive, Black, Sterile, Blister	
		Packed	
	332112	200 µl Automation Filter Tips for	96 Tips*2 Blisters/pk, 2304/cs
		Tecan, Conductive, Black, Sterile,	
		Blister Packed	
and the second	332202	1000 µl Automation Tips for Tecan,	96 Tips*2 Blisters/pk, 2304/cs
		Conductive, Black, Sterile, Blister	
		Packed	
	332212	1000 µI Automation Filter Tips for	96 Tips*2 Blisters/pk, 2304/cs
		Tecan, Conductive, Black, Sterile,	
		Blister Packed	
	332016	50 µl Automation Filter Tips for Tecan,	96 /pk, 4800/cs
	332010	Conductive, Black, Sterile, Box-packed	
	332116	200 µl Automation Filter Tips for	96 /pk, 4800/cs
	332110	Tecan, Conductive, Black, Sterile,	

 Head Office
 No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China

 Tel:
 +86+ 510-6800 6788
 Email: info@nest-wuxi.com

 Online:
 www.cell-nest.com
 Email: info@nest-wuxi.com

Oversea







	Box-packed	
	1000 µI Automation Filter Tips for	96 /pk, 4800/cs
332216	Tecan, Conductive, Black, Sterile,	
	Box-packed	

	A	Automation Tips for Tecan, Clear	
	332003	50 μl Automation Tips for Tecan, Clear, Sterile, Blister Packed	96 Tips*2 Blister/pk, 2304/cs
	332013	50 μl Automation Filter Tips for Tecan, Clear, Sterile, Blister Packed,	96 Tips*2 Blisters/pk, 2304/cs
And a state of the	332103	200 µl Automation Tips for Tecan, Clear, Sterile, Blister Packed	96 Tips*2 Blisters/pk, 2304/cs
	332113	200 μl Automation Filter Tips for Tecan, Clear, Sterile, Blister Packed	96 Tips*2 Blisters/pk, 2304/cs
	332203	1000 μl Automation Tips for Tecan, Clear, Sterile, Blister Packed,	96 Tips*2 Blisters/pk, 2304/cs
	332213	1000 μl Automation Filter Tips for Tecan, Clear, Sterile, Blister Packed	96 Tips*2 Blisters/pk, 2304/cs
	332007	50 μl Automation Tips for Tecan, Clear, Sterile, Racked,	96/pk, 4800/cs
	332017	50 μl Automation Filter Tips for Tecan, Clear, Sterile, Racked	96/pk, 4800/cs
	332107	200 μl Automation Tips for Tecan, Clear, Sterile, Racked	96/pk, 4800/cs
	332117	200 µl Automation Filter Tips for Tecan, Clear, Sterile, Racked	96/pk, 4800/cs
	332207	1000 μl Automation Tips for Tecan, Clear, Sterile, Racked	96/pk, 4800/cs
	332217	1000 μl Automation Filter Tips for Tecan, Clear, Sterile, Racked	96/pk, 4800/cs

Automation Conductive Tips for Hamilton, Black			
	245002	50 μl Automation Tips for Hamilton,	Blister Packed, 96/pk, 2304/cs
	345002	Conductive, Black, Sterile	
	345012	50 µl Automation Filter Tips for	Blister Packed, 96/pk, 2304/cs
		Hamilton, Conductive, Black, Sterile	
	345102	300 µl Automation Tips for Hamilton,	Blister Packed, 96/pk, 2304/cs
		Conductive, Black, Sterile	
	345112	300 µl Automation Filter Tips for	Blister Packed, 96/pk, 2304/cs
Manual Maria			

Head Office

No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com Online: www.cell-nest.com

Oversea







	1	1505001 6	(13013485 & 13011137
		Hamilton, Conductive, Black, Sterile	
	345202	1000 µl Automation Tips for Hamilton,	Blister Packed, 96/pk, 2304/cs
		Conductive, Black, Sterile	
	245212	1000 µl Automation Filter Tiprs for	Blister Packed, 96/pk, 2304/cs
	345212	Hamilton, Conductive, Black, Sterile	
		50 µl Automation Filter Tips for	Blister Packed, 96/pk, 2304/cs
L HINK CHANNE	345013	Hamilton, Conductive, Black, with	
		Barcode, Sterile	
		300 µl Automation Filter Tips for	Blister Packed, 96/pk, 2304/cs
	345113	Hamilton, Conductive, Black, with	
A CONTRACT OF CONTRACT OF		Barcode, Sterile	
		1000 µl Automation Filter Tiprs for	Blister Packed, 96/pk, 2304/cs
	345213	Hamilton, Conductive, Black, with	
		Barcode, Sterile	
		50 µl Automation Filter Tips for	Blister Packed, 96 Tips*5
	345015	Hamilton, Conductive, Black, with	Blisters/pk, 5760/cs
and the second se		Barcode, Sterile,	
	345005	50 µl Robotic Tips for Hamilton,	96 Tips*5 Blisters/pk, 5760/cs
	345005	Conductive, with Barcode, Sterile	
		300 µI Automation Filter Tips for	96 Tips*5 Blisters/pk, 5760/cs
	345115	Hamilton, Conductive, Black, with	
		Barcode, Sterile, Blister Packed	
	345105	300 μl Robotic Tips for Hamilton,	96 Tips*5 Blisters/pk, 5760/cs
	345105	Conductive, with Barcode, Sterile	
		1000 µI Automation Filter Tips for	Blister Packed, 96 Tips*5
	345215	Hamilton, Conductive, Black, with	Blisters/pk, 3840/cs
		Barcode, Sterile,	
	045005	1000 µl Robotic Tips for Hamilton,	96 Tips*5 Blisters/pk, 3840/cs
	345205	Conductive, with Barcode, Sterile	

Automation Tips for Hamilton, Clear				
	345503	50 µl Automation Tips for Hamilton,	Blister Packed, 96/pk, 2304/cs	
		Clear, Sterile, with Barcode		
	345513	50 µl Automation Filter Tips for	Blister Packed, 96/pk, 2304/cs	
A CONTRACT OF		Hamilton, Clear, Sterile, with Barcode		
Martin Standy of the	345603	300 μl Automation Tips for Hamilton,	Blister Packed, 96/pk, 2304/cs	
		Clear, Sterile, with Barcode		
	345613	300 µI Automation Filter Tips for	Blister Packed, 96/pk, 2304/cs	
		Hamilton, Clear, Sterile, with Barcode,		
	345703	1000 µl Automation Tips for Hamilton,	Blister Packed, 96/pk, 2304/cs	

 Head Office
 No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China

 Tel:
 +86+ 510-6800 6788
 Email: info@nest-wuxi.com

 Online:
 www.cell-nest.com
 Mail: info@nest-wuxi.com

Oversea

DATA SHEET






ISO9001 & ISO13485 & ISO11137

		Clear, Sterile, with Barcode	
	345713	1000 µl Automation Filter Tips for	Blister Packed, 96/pk, 2304/cs
		Hamilton, Clear, Sterile, with Barcode	

Automation Tips for Beckman, Clear			
	317511	250 μl Automation Tips for Beckman, Clear, Sterile	Racked, 96/pk, 4800/cs
	317501	190 μl Automation Filter Tips for Beckman, Clear, Sterile	Racked, 96/pk, 4800/cs
	317011	1070 µl Automation Tips for Beckman, Clear, Sterile	Racked, 96/pk, 4800/cs
	317101	1000 µl Automation Filter Tips for Beckman, Clear, Sterile	Racked, 96/pk, 4800/cs

Technical Drawing

T-brand Automation Tips



 Head Office
 No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
 Oversea


 Tel:
 +86+ 510-6800 6788
 Email: info@nest-wuxi.com

 Online:
 www.cell-nest.com




CE

H-brand Automation Tip (Filter)




B-brand Automation Tips



Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com Online: www.cell-nest.com





Tecan Conductive Tip Rack

 Head Office
 No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
 Oversea
 NEST USA (New Jersey/ Phoenix) NEST scientific 株式会社 (Yokohama,Japan)

 Tel: +86+ 510-6800 6788
 Email: info@nest-wuxi.com
 NEST Scientific Europe B.V (Netherlands)

 Online: www.cell-nest.com
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)